

CHEMICAL THREAT PACKAGING AND SHIPPING GUIDE FOR BLOOD SPECIMENS

COLLECTED FROM INDIVIDUALS WHO MAY HAVE BEEN EXPOSED TO CHEMICAL AGENTS

II. PACKAGING SPECIMENS – BLOOD TUBES

INNER PACKAGING	OUTER PACKAGING
<p>Package all blood tubes from the same patient together.</p> <ol style="list-style-type: none"> 1. Place absorbent material between the blood tubes and the 1st layer of secondary packaging. 2. Use enough absorbent material to absorb the entire contents of the blood tubes. 3. Separate or wrap tubes to prevent tube-to-tube contact. 4. Secure the 1st layer of secondary packaging with one continuous strip of tamper-evident tape and write packer's initials half on the tape and half on the 1st layer of secondary packaging. <ol style="list-style-type: none"> a. Pack blood tubes in a gridded box lined with absorbent material—or—pack a sealable polystyrene foam container or blood tube shipment sleeve and transport tube with individually wrapped tubes. b. Seal the box or polystyrene foam container or transport tube with one continuous piece of evidence tape and write your initials half on the tape and half on the container. 5. Wrap and seal the first layer of secondary packaging (e.g., gridded box) with absorbent material. 6. Seal one wrapped gridded box or alternative container inside a clear, leak-proof biohazard polybag equivalent to Saf-T-Pak product STP-701, STP-711 or STP-731. 7. Place this bag inside a white Tyvek® outer envelope (or equivalent); seal the opening with a continuous strip of evidence tape initialed half on the packaging and half on the tape. 8. According to 49 CFR 173.199(b), for air transport, either the primary receptacle or the secondary packaging must be capable of withstanding, without leaking, an internal pressure producing a pressure differential of not < 95 kPa (0.95 bar, 14 psi). 	<ol style="list-style-type: none"> 1. Use a polystyrene foam-insulated, corrugated fiberboard shipper (may be available from your transfusion service or send-outs department). 2. For cushioning, place additional absorbent material in the <u>bottom</u> of the shipper. 3. Add a single layer of refrigerator packs on top of absorbent material. 4. Place the packaged specimens on top of the refrigerator packs. 5. Use additional cushioning material to minimize shifting while the shipper is in transit. 6. Place additional refrigerator packs on top of the secondary packaging to maintain a shipping temperature of 1° C - 10° C for the duration of transit. 7. Place blood shipping manifest in a sealable plastic bag and put on top of packs inside the shipper. 8. Keep chain-of-custody documents for your files. 9. Place lid on shipper and secure corrugated fiber board (box) with filamentous shipping tape. 10. Place your return address (including telephone number) in the upper left-hand corner of the shipper top and put MDPH's receiving address in center. 11. Affix labels and markings adjacent to the shipper's/consignee's address that appears on the shipper. 12. Place the UN 3373 label and the words "Biological Substance, Category B" adjacent to the label on the front of the shipper. 13. SHIP TO: William A. Hinton State Laboratory Institute Chemical Threat Response Laboratory 305 South Street, Room 305 Jamaica Plain, MA 02130 ATTN: Dr. Jennifer Jenner 24/7 Cell: 617-839-1283 Tel: 617-983-6650

Massachusetts Department of Public Health
William A. Hinton State Laboratory Institute
Chemical Threat Response Laboratory
305 South Street, Jamaica Plain, MA 02130
Tel: 617-983-6650 Fax: 617-983-6662