ATTACHMENT 11



The Commonwealth of Massachusetts

Executive Office of Health and Human Services
Department of Public Health
State Laboratory Institute
305 South Street, Jamaica Plain, MA 02130

MITT ROMNEY GOVERNOR

KERRY HEALEY LIEUTENANT GOVERNOR

RONALD PRESTON SECRETARY

CHRISTINE C. FERGUSON COMMISSIONER

To: Local Boards of Health

Animal Control Officers and Inspectors

From: Barbara Werner, Ph.D., Director

Infectious Disease Laboratories Division, Massachusetts Department of Public Health

Date: March 2003

Re: Revised Guidelines for Specimen Transport and Rabies Testing

The Massachusetts Department of Public Health (MDPH) Bureau of Laboratory Sciences is committed to providing accurate, timely testing of suspect rabid animals involved in exposures of humans or domestic animals. Together with the Bureau of Communicable Disease Control (BCDC), we will continue to provide quality rabies testing services and information for Commonwealth residents.

The U.S. Department of Transportation (DOT) recently enhanced its packaging requirements for diagnostic and infectious specimens. Please see the attached Guidelines for Specimen Transport and Rabies Testing, which replace previous versions of these Guidelines. The revised Specimen Request Form for Rabies Testing is enclosed. Please begin using this form immediately. The Guidelines and Specimen Request Form for Rabies Testing are also available on the MDPH website at http://www.state.ma.us/dph/bls/viro.htm.

For specific enhancements to DOT regulations concerning the transportation of diagnostic specimens, please consult the Federal Register, Vol. 67, No. 157. Part III, Department of Transportation, Research and Special Programs, 49 CFR Part 171 et al., Hazardous Materials: Revision to Standards for Infectious Substances; Final Rule, available at http://hazmat.dot.gov/67fr-53118.pdf.

Questions concerning testing and specimen transport should be directed to the Virology Laboratory at (617) 983-6385, Monday through Friday between 8am and 6pm. Questions concerning exposures, the need for testing, and prophylaxis should be directed to the Division of Epidemiology and Immunization, BCDC, MDPH, at (617) 983-6800.

Specimen Request Form for Rabies Testing STATE LABORATORY INSTITUTE

305 South Street

Jamaica Plain, MA 02130-3597

Do not use this space								
Date Received	1	1						

PLEASE PRINT		Tel. 617-	983-6385	DO NOT				
1. PROVIDER/SENDER INFORMATION			2. OWNER INFORMATION (or person who found animal)					
Name				Name: Last		Firs	t Initial	
Address: No./Street/Apt.#	‡			Address: No./	/Street/Apt #			
City/Town	State Zip Code		City/Town State Zip Code					
Phone Number: ()			Phone Number	er: ()				
3. SPECIMEN INFOR	MATION	□ Pet	□ Stray	□ Wild	□ Unknov	wn		
Species	Bre	ed	Age	/			Natural Euthanized Method	
Location Where	Symptom	ıs:	□ Found D			Reason for Rabies Testing:		
Animal was Found:	□ Aggr		Lethargy		Jnexplained			
Town	□ Atax		Paralysis		Wound		xposure	
	– 🛚 Diso	rientation	Salivation	n 🗆 (Other	Acting	g Sick	
County_		D:44	L A 41	. A ! I ! D	1 40		18-4	
Travel Out of State:			en by Another	Animai in Pa	ast 12	Vaccination History: date □ Rabies Vaccinated (/ /)		
☐ Yes (Location_	Data /		nths:	nimal	,			
Location	Date/	// 🗀	No	animal) □ Not Rabies Vaccinated (not current) □ Unknown				
☐ Unknown			Unknown			- CHRIO	vvii	
4. EXPOSURE INFORMATION								
Person(s) Exposed Exposure Date / /			Animal(s) E	Exposed	Ex	posure Date//		
Name Physician Name			Name					
Address: No./Street/Apt.#			Species Age					
City/Town State Zip Code			Address: No./Street/Apt.# (if different from owner)					
Phone Number: ()	Physic	ian Number: ()	City/Town		State	Zip Code	
Type of □	Bite	Body	Site			Bite	Body Site	
	Scratch			T	П	Scratch	-	
	Lick			Type of Exp (check one)		Lick		
		Seve	rity	(Check one)	_	Other	Severity	
	Unknown			***************************************		Unknown		
Circumstance of Exposure:	Capture		Specimen Preparation	0:		Fight		
of Exposure:	Unprovoked Att Provoked Attack		Other	Circumstar Exposure:	ice of \Box	Vicinity	nal Contact	
	Handling			(check one)		Other	nai Contact	
				(Cricck one)		Otrici		
Pre-Exposure Vaccir			sure Rabies	Pre-Exposu	ire Vaccinatio	n History:	Post-Exposure Rabies	
History:		Vaccinatio	n:	□ Rabies Vaccinated (current) Vaccination:			Vaccination:	
□ Rabies Vaccina	ted (current)	□ Rece	ived	□ Not Rabies Vaccinated □ Received				
□ Not Rabies Vac			Received	(not current)				
(not current)		□ Unkn		☐ Unknown ☐ Unknown				
□ Ùnknown ´								
5. FLUORESCENT R	ABIES ANTIB	ODY TEST	RESULTS: (for	lab use only)	Reported By	(initials)	Date://	
□ Positive (rabid) □ Negative (not rabid) □ Specimen								
Unsatisfactory Comments:								

Guidelines for Sending Specimens for Rabies Testing

Recently the U.S. Department of Transportation enhanced its packaging requirements for specimens. Federal regulations require that diagnostic specimens be triple packaged according to specific guidelines for transport. The guidelines below are for shipping specimens for rabies testing to the State Laboratory Institute (SLI). A Specimen Request Form for Rabies Testing must be included with each specimen.

Materials:

- 2 zipper-lock or otherwise sealable, leak-proof plastic bags
- absorbent material (pads, paper towels, etc.)
- 1 frozen cool pack
- 1 pair disposable gloves

- 1 permanent marker
- 1 cardboard box or plastic foam mailing cooler
- cushioning material (newspaper, etc.)

Specimen Packaging and Transport:

- 1. All animals involved in exposures of humans or domestic animals will be tested. With the exception of bats, the head **only** should be prepared for transport. Please note that "road kills" may not be appropriate for testing and should not be sent if the head is severely damaged.
- 2. With gloved hands, place the specimen into a plastic bag. Seal the bag.
- 3. Place the bagged specimen into the larger plastic bag with enough absorbent material to absorb any leaking fluid.
- 4. Place a frozen cool pack into the bag. DO NOT USE ICE. Seal the bag.
- 5. Place the wrapped, double-bagged specimen inside a cardboard box or cooler with adequate cushioning material. Place the completed Specimen Request Form for Rabies Testing in the cardboard box.
- 6. Seal the box. Label the box "Diagnostic Specimen" in permanent marker. Clearly label with the destination address (State Laboratory Institute, Rabies Laboratory, Room 710, 305 South St., Jamaica Plain, MA, 02130), and the sender's address in permanent marker.

Specimen Transport:

Specimens may be delivered to SLI via private courier, U.S. Postal Service, United Parcel Service (UPS), and FedEx. FedEx will not accept whole animals for shipment, and requires that packages be at least 7" x 4" x 2" in size. Please contact the shipping agency for questions regarding its requirements for packaging. DO NOT leave specimens in delivery service drop boxes.

Specimens should be delivered to SLI as soon as possible. The specimen must be maintained cold, but not frozen at all times. Specimens received more than 24 hours after death or after freezing may not be suitable for testing.

Testing and Exposure:

Animals involved in a human or domestic animal exposure will be tested as soon as possible. Results will be phoned to the sender as soon as the test is completed. Specimens received by 1:00 pm will be prepared and tested the same day. The Bureau of Communicable Disease Control (BCDC) epidemiologists can authorize emergency testing on weekends.

Questions concerning testing and specimen transport should be directed to the Virology Laboratory at (617) 983-6385, Monday through Friday between 8am and 6pm. Questions concerning exposures, need for testing, and prophylaxis should be directed to the Division of Epidemiology and Immunization, BCDC, MDPH, at (617) 983-6800.