



Title: Quality Assurance and Quality Control Requirements and Performance Standards for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

WSC – CAM – III B

Quality Assurance and Quality Control Requirements for ***SW-846 Methods 7470A and 7471A, Mercury in Liquid Waste and Mercury in Solid or Semi-Solid Waste (Manual Cold-Vapor Technique)*** for the Massachusetts Contingency Plan (MCP)

Document Status: Final

Donald Muldoon
Department of Environmental Protection
One Winter Street
Boston, MA 02108
donald.muldoon@state.ma.us

Document Chronology

Revision	Status	Date
0	Public Comment Draft	19 January 01
1	Revised Public Comment Draft	26 June 01
2	Interim Final	31 August 01
3	Revised Interim Final	12 July 2002
4	Final	1 May 03



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

III. Metal Analysis

B. Quality Assurance and Quality Control Requirements and Performance Standards for SW-846 Methods 7470A and 7471A, Mercury in Liquid Waste and Mercury in Solid or Semi-Solid Waste (Manual Cold-Vapor Technique)

B.1 TABLE OF CONTENTS

1.0	QA/QC Requirements for SW-846 Methods 7470A and 7471A	3
1.1	Overview of Methods	3
1.2	General QA/QC Requirements for SW-846 Methods 7470A and 7471A	5
1.3	Specific QA/QC Requirements and Performance Standards for SW-846 Methods 7470A and 7471A	6
2.0	Data Usability Assessment for SW-846 Methods 7470A and 7471A	12
3.0	Reporting Requirements for SW-846 Methods 7470A and 7471A	12
3.1	General Reporting Requirements for SW-846 Methods 7470A and 7471A	12
3.2	Specific Reporting Requirements for SW-846 Methods 7470A and 7471A	12
4.0	Regulatory Limits for Mercury under 310 CMR 40.000	13

List of Tables, and Appendices

Number	Title	Page
Table III B-1	Performance Requirements for SW-846 Methods 7470A and 7471A	8-11
Table III B-2	Routine Analytical Reporting Requirements for SW-846 Methods 7470A and 7471A	13
Appendix III B-1	Sample Collection, Preservation and Handling Procedures for Mercury in Liquid Waste and Mercury in Solid or Semi-Solid Waste (Manual Cold-Vapor Technique)	14
Appendix III B-2	Minimum Field QC Sample Frequency in Support of "Presumptive Certainty" Status for Analyzed for Metals by CVAA (SW-846 methods 7470A/7471A)	15



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

1.0 QA/QC Requirements for SW-846 Methods 7470A and 7471

1.1 Overview of Methods

1.1.1 SW-846 Method 7470A MERCURY IN LIQUID WASTE

SW-846 Method 7470A is a cold-vapor atomic absorption (CVAA) procedure for determining the concentration of mercury in mobility-procedure extracts, aqueous wastes, surface and ground waters. Prior to analysis, the liquid samples must be pretreated according to the procedure described in Section 7.0 of this method.

Quantitation of Mercury is based on the absorption of radiation at 253.7-nm by mercury vapor. The mercury is reduced to the elemental state and aerated from solution in a closed system. The mercury vapor passes through a cell positioned in the light path of an atomic absorption spectrophotometer. Absorbance (peak height) is measured as a function of mercury concentration. The typical aqueous detection limit for this method is 0.0002 mg/L (0.2 µg/L)

Sample preservation, container and analytical holding time specifications for surface water, groundwater, soil, and sediment matrices for metals analyzed in support of MCP decision-making are presented in Appendix III B–1 of this document and Appendix VII-A, CAM–VII A, “Quality Assurance and Quality Control Guidelines for the Acquisition and Reporting of Analytical Data in Support of Response Actions Conducted Under the Massachusetts Contingency Plan (MCP)”.

1.1.2 SW-846 Method 7471A MERCURY IN SOLID OR SEMISOLID WASTE

SW-846 Method 7471A is also a CVAA procedure based on the absorption of radiation at 253.7-nm wavelength by mercury vapor, approved for determining total mercury (organic and inorganic) in soils, sediments, bottom deposits, and sludge-type materials. For routine MCP mercury analyses in solid or semisolid waste, a 0.2-gram sample (dry-weight basis) is subjected to an appropriate dissolution step, as described in Section 7.0 of the Method, prior to analysis. It should be noted, that if the selected dissolution procedure is not sufficient to dissolve a specific matrix type or sample, then this method is not applicable for that matrix.

Because of this small sample (0.2-grams, as recommended in Section 7.1, Procedure), “as-received” field samples of solid or semisolid materials must be thoroughly homogenized in the laboratory prior to mercury analysis to obtain a representative aliquot for analysis. If mercury is a contaminant of concern (COC) at an MCP site or if poor precision and/or accuracy (that would adversely effect MCP decision making) associated with sample matrix heterogeneity is anticipated, then the following preemptive or corrective measures should be considered by the LSP:



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

- Direct the laboratory to use a larger sample aliquot (up to 10 grams) to enhance sensitivity and/or precision. The laboratory must adjust the concentration/volume of all reagents to accommodate the analytical requirements of the larger sample aliquot and re-demonstrate applicable laboratory-specific QC control limits and performance criteria;
- Direct the laboratory to prepare and analyze replicate 0.2-gram aliquots (two or more) for each field sample to better assess variability, associated with sample matrix heterogeneity;
- Utilize more effective field sample homogenization procedures (prior to submitting sample to laboratory for mercury analysis); and/or
- Increase frequency of field duplicates.

If overall site heterogeneity for mercury is a specific concern (i.e. distribution of mercury contamination is inconsistent with the conceptual site model), it is recommended that the number of field samples be increased to improve the representativeness of mercury results for site assessment and characterization.

SW-846 Method 7471A Mercury In Solid Or Semisolid Waste Analytical Notes:

For routine MCP Mercury analysis of solid or semisolid samples, the Department does not require the analysis of triplicate sample aliquots for each determination as described in SW-846 Method 7471A, Section 7.1, Sample Preparation. Consistent with the requirements of 310 CMR 40.0017(1), mercury analysis must be scientifically valid and defensible, and of a level of precision and accuracy commensurate with its stated or intended use. The requirement for triplicate analysis for every determination is overly burdensome for most MCP applications. However, replicate analyses may be warranted in some cases to evaluate sample matrix heterogeneity, at the discretion of the LSP.

1.1.3 SW-846 Methods 7470A and 7471A Interferences

Samples submitted to a laboratory for trace metal analysis may become contaminated by numerous routes during both sampling and analysis. Potential sources of contamination may include:

- Metallic or metal-containing containers and sampling equipment,
- Laboratory acids or reagents,
- Improperly cleaned or stored equipment , and
- Atmospheric inputs such as dirt and dust

Other potential contaminants and their recommended corrective measures, as appropriate, are listed below:



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

- Potassium permanganate is added to eliminate possible interference from sulfide. Concentrations as high as 20 mg/Kg of sulfide, as sodium sulfide do not interfere with the recovery of inorganic mercury added to reagent water.
- Seawaters, brines, and industrial effluents high in chlorides require additional permanganate (as much as 25 mL) because, during the oxidation step, chlorides are converted to free chlorine, which also absorbs radiation of 253.7 nm. Care must therefore be taken to ensure that free chlorine is absent before the mercury is reduced and swept into the cell. This may be accomplished by using an excess of hydroxylamine sulfate reagent (25 mL).
- Certain volatile organic materials that absorb at the optimum analytical wavelength (253.7 nm) may also cause interference at elevated concentrations. A preliminary run without reagents should determine if this type of interference is present.
- Copper has also been reported to interfere; however, copper concentrations as high as 10 mg/Kg have been shown to have no effect on recovery of mercury from spiked samples.

1.2 General QA/QC Requirements for SW-846 Methods 7470A and 7471A

Each laboratory that uses SW-846 Methods 7470A and 7471A is required to operate a formal quality assurance program to demonstrate the precision and bias of the method as performed by the laboratory and procedures for determining the method reporting limit (RL). The minimum requirements of this program consist of an initial demonstration of laboratory proficiency, ongoing analysis of standards and blanks as a test of continued performance, and the analysis of laboratory control spikes (LCSs) and matrix spikes (MS), to assess accuracy and/or precision. Matrix duplicates or matrix spike duplicates (MSD) may also be used to evaluate precision when such samples are analyzed either at discretion of laboratory or at request of data-user.

Laboratories must document and have on file an Initial Demonstration of Proficiency for each combination of sample preparation and determinative method being used. These data must meet or exceed the performance standards as presented in Section 1.3 and Table III B-1 of this method. Procedural requirements for performing the Initial Demonstration of Proficiency can be found in SW-846 Chapter One, Section 4.4.1 and SW-846 Methods 7470A and 7471A, Section 8.0. The data associated with the Initial Demonstration of Proficiency must be kept on file at the laboratory and made available to potential data-users on request. The data associated with the Initial Demonstration of Proficiency for SW-846 Method 7470A and 7471A must include the following:



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

QC Element	Performance Criteria
Initial Calibration	CAM-III B, Table III B-1
Continuing Calibration	CAM-III B, Table III B-1
Method Blanks	CAM-III B, Table III B-1
Percent Recovery for MS/LCS	CAM-III B, Table III B-1
Percent Relative Standard Deviation for Duplicates	SW-846 Method Chapter One, Section 4.4.3

It is essential that laboratory-specific performance criteria for LCS and the other data quality indicators, listed in Table III B-1, also be calculated and documented. When experience indicates that the criteria recommended in specific methods are frequently not met for some analytes and/or matrices, the in-house performance criteria will be a means of documenting these repeated exceedances. Laboratories are encouraged to actively monitor pertinent quality control performance standards described in Table III B-1 to assess analytical trends (i.e., systematic bias, etc) and to improve overall method performance.

For SW-846 Methods 7470A and 7471A, laboratory-specific control limits must meet or exceed (demonstrate less variability than) the performance standards for each QC element listed on Table III B-1. It should be noted that the performance standards listed in Table III B-1 are based on multiple-laboratory data, which are in most cases expected to demonstrate more variability than performance standards developed by a single laboratory. Laboratories are encouraged to continually strive to minimize variability and improve the accuracy and precision of their analytical results. The LSP will utilize this analytical performance data to verify that the results reported by the laboratory are consistent with the pre-established data quality objectives for the disposal site.

Use of this method is restricted to use by, or under the supervision of, analysts who are knowledgeable of cold-vapor atomic absorption spectrophotometry as a quantitative tool and the correction of spectral, chemical, and physical interferences described in this method.

1.3 Specific QA/QC Requirements and Performance Standards for SW-846 Methods 7470A and 7471A

Specific QA/QC requirements and performance standards for SW-846 Methods 7470A and 7471A are presented in Table III B-1. Strict compliance with the QA/QC requirements and performance standards for this method, as well as satisfying other analytical and reporting requirements will provide an LSP with "Presumptive Certainty" regarding the usability of analytical data to support MCP decisions. The concept of "Presumptive Certainty" is explained in detail in Section 2.0 of CAM VII-A.



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

While optional, parties electing to utilize these protocols will be assured of “Presumptive Certainty” of data acceptance by agency reviewers. In order to achieve “Presumptive Certainty”, parties must:

- Comply with the procedures described and referenced in CAM–III B;
- Comply with the applicable QC requirements prescribed in Table III B-1 for this test procedure;
- Collect and analyze field QC samples at the frequencies prescribed in the CAM–VII A and Appendix III B-2 for this test procedure; and
- Adopt the reporting formats and elements specified in the CAM–VII A.

In achieving the status of “Presumptive Certainty”, parties will be assured that analytical data sets:

- ✓ Will satisfy the broad QA/QC requirements of 310 CMR 40.0017 and 40.0191 regarding the scientific defensibility, precision and accuracy, and reporting of analytical data;
- ✓ May be used in a data usability assessment, and if in compliance with all MCP Analytical Method standards, laboratory QC requirements, and field QC recommended limits and action levels, the data set will be considered useable data to support site characterization decisions made pursuant to the MCP; and
- ✓ May be used to help support a data representativeness assessment.

Widespread adherence to the “Presumptive Certainty” approach will promote inter-laboratory consistency and provide the regulated community with a greater degree of certainty regarding the quality of data used for MCP decision-making. The issuance of these requirements and standards is in no way intended to preempt the exercise of professional judgement by the LSP in the selection of alternative analytical methods. However, parties who elect not to utilize the “Presumptive Certainty” option have an obligation, pursuant to 310 CMR 40.0017 and 40.0191(2)(c), to demonstrate and document an overall level of (laboratory and field) QA/QC, data usability, and data representativeness that is adequate for and consistent with the intended use of the data



Title: **Specific QA/QC Requirements and Performance Standards for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste**

Required QA/QC	Data Quality Objective	Performance Standard	Required Deliverable	Recommended Corrective Action	Analytical Response Action
Preparation of Samples	Accuracy and Representativeness	All aqueous and solid samples must be prepared (digested) prior to analysis. See Methods 7470 and 7471 for details. Note: MADEP requires only 1 preparation for each field sample to generate a reportable result.	No	See section 2.1.2 in this guidance for further details on obtaining representative soil results.	
Initial Calibration	Laboratory Analytical Accuracy	(1) Frequency - Daily prior to sample analysis (2) Minimum of a calibration blank plus five calibration standards (3) Linear curve fit with correlation coefficient $r \geq 0.995$. Second order curve fit may be used if $r \geq 0.995$.	No	Re-calibrate as required by method.	
Initial Calibration Verification (ICV)	Laboratory Analytical Accuracy	(1) Frequency - Daily immediately after initial calibration (2) Separate-source from calibration standards (3) ICV % recovery must be 90-110%;	No	Re-calibrate/Re-analyze ICV as required by method	Suspend all analyses until Initial Calibration non-conformance is rectified.
Initial Calibration Blank (ICB)	Laboratory Analytical Sensitivity (instrument drift and contamination evaluation)	(1) Frequency - Daily immediately after ICB (2) Must be matrix-matched (the same concentration of acids as standards and samples) (3) ICB must be < RL	No	Re-calibrate/Re-analyze ICB as required by method	



Title: **Specific QA/QC Requirements and Performance Standards for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste**

Required QA/QC	Data Quality Objective	Performance Standard	Required Deliverable	Recommended Corrective Action	Analytical Response Action
Low Level Calibration Check Standard	Laboratory Analytical Sensitivity (to support the RL)	Only required if NOT including the RL as the low-level standard in the calibration curve. (1) Frequency - Daily prior to field sample analysis (2) Check Standard at the level of the RL for mercury (3) % recovery recommended 70-130%	No	Re-calibrate / narrate	Narrate non-compliance.
Continuing Calibration Verification (CCV)	Laboratory Analytical Accuracy	(1) Frequency - Every 10 samples and at end of run (2) Same-source as calibration standards; near mid-point of linear range (3) CCV % recovery must be 80-120%	No	Re-calibrate/Re-analyze all samples since last compliant CCV.	Narrate non-compliance.
Continuing Calibration Blank (CCB)	Laboratory Analytical Sensitivity (instrument drift and contamination evaluation)	(1) Frequency - Every 10 samples and at end of run immediately after CCV (2) Must be matrix-matched (the same concentration of acids as standards and samples) (3) CCB must be < RL	No	Re-calibrate/Re-analyze all samples since last compliant CCB.	Narrate non-compliance.
Method (Preparation) Blank	Laboratory Method Sensitivity (contamination evaluation)	(1) Frequency - One per digestion batch of < 20 field samples. (2) Must be matrix-matched (the same concentration of acids as calibration and QC standards) and digested with the samples (3) Method Blank must be < RL	Yes	Re-digest/Re-analyze all associated samples unless all detected results are > 10x method blank level.	Narrate non-compliance.



Title: **Specific QA/QC Requirements and Performance Standards for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste**

Required QA/QC	Data Quality Objective	Performance Standard	Required Deliverable	Recommended Corrective Action	Analytical Response Action
Laboratory Control Sample (LCS)	Laboratory Method Accuracy	(1) Frequency - One per digestion batch of < 20 field samples. (2) LCS must be matrix-matched (aqueous or solid) to field samples, a separate-source from the calibration standards, and digested with the samples. (3) LCS % recovery for mercury must be 80-120% for aqueous media and vendor control limits (95% confidence limits) for solids.	Yes	Re-digest/Re-analyze all associated samples.	Narrate non-compliance. Note: MADEP modification to frequency of LCS for consistency with other methods.
Matrix Spike Sample (MS)	Method Accuracy in Sample Matrix	(1) Frequency - One per digestion batch of < 20 field samples. (2) MS % recoveries must be between 75 – 125 for all matrices (3) Laboratories are expected to develop their own in-house control limits for each media, which should fall within the limits listed above.	Yes	No corrective action required.	Narrate non-compliance
Matrix Spike Duplicate Sample (MSD)	Method Precision in Sample Matrix	(1) Frequency - One per digestion batch of < 20 field samples. (2) MSD relative percent difference (RPD) criteria: aqueous results: $\pm 20\%$; soil and sediment results: $\pm 35\%$.	Yes	No corrective action required.	Narrate non-compliance



Title: **Specific QA/QC Requirements and Performance Standards for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste**

Required QA/QC	Data Quality Objective	Performance Standard	Required Deliverable	Recommended Corrective Action	Analytical Response Action
Matrix Duplicate Sample (MD)	Method Precision in Sample Matrix	(1) Optional - One per digestion batch of < 20 field samples. (2) MD relative percent difference (RPD) criteria recommended: aqueous results > 5x RL: $\pm 20\%$; aqueous results < 5x RL: difference \leq RL; soil and sediment results > 5x RL: $\pm 35\%$; soil and sediment results < 5x RL: difference $\leq 2x$ RL.	Yes, if MD requested by LSP or run by laboratory as routine QA/QC	No corrective action required,	Narrate non-compliance. Note: poor sample/MD precision may be an indication of sample matrix heterogeneity. See Section 2.1.2 of this guidance for recommended procedures to obtain representative mercury data in soil/sediment matrices.
General Reporting	NA	(1) Non detected values must be reported with the sample-specific reporting limit for each ICP analyte. (2) The RL must be \leq the applicable regulatory compliance standard for each metal reported (3) The RL must be verified at least daily with a low-level calibration check standard following the calibration curve or supported by the low-level standard in the calibration curve. (4) Results for soils/sediments must be reported on a dry weight basis for comparison to MCP regulatory standards (5) Sample concentrations that exceed the highest calibration standard must be diluted (in the same acid matrix) to fall within the calibration range when reanalyzed	Yes	Not applicable.	



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

2.0 Data Usability Assessment for SW-846 Methods 7470A/7471A

Overall data usability is influenced by uncertainties associated with both sampling and analytical activities. This document provides detailed quality control requirements and performance standards for SW-846 Methods 7470A/7471A which may be used to assess the analytical component of data usability. The sampling component of data usability, an independent assessment of the effectiveness of sampling activities to meet data quality objectives, is not substantively addressed in this document.

Note to data users: Consistent with USEPA Region I data validation guidance, the Department requires rejection of non-detected metals results with <30% MS recovery. If Mercury MS % recovery is < 30%, and non-detected results were found, the LSP should consider the use of one or more of the following actions to assess and /or overcome the sample matrix effect encountered:

- 1) Re-homogenize sample. Prepare and analyze a new sample/MS pair;
- 2) Perform method of standard additions for quantitation;
- 3) Perform serial dilution; and/or
- 4) Perform post-digestion spike.

3.0 Reporting Requirements for SW-846 Methods 7470A and 7471A

3.1 General Reporting Requirements for SW-846 Methods 7470A and 7471A

General reporting requirements for analytical data used in support of assessment and evaluation decisions at MCP disposal sites are presented in WSC-CAM-VIIA. This guidance document provides recommendations for field QC, as well as the required content of the Environmental Laboratory Report, including

- Laboratory identification information presented in CAM-VII A, Section 2.4.1,
- Analytical results and supporting information in CAM-VII A, Section 2.4.2,
- Sample- and batch-specific QC information in CAM-VII A, Section 2.4.3,
- Laboratory Report Certification Statement in CAM-VII A, Section 2.4.4,
- Copy of the Analytical Report Certification Form in CAM-VII A, Exhibit VII A-1,
- Environmental Laboratory Case Narrative contents in CAM-VII A, Section 2.4.5,
- Chain of Custody Form requirements in CAM-VII A, Section 2.4.6

3.2 Specific Reporting Requirements for SW-846 Methods 7470A and 7471A

Specific analytical reporting requirements for SW-846 Methods 7470A/7471A are presented in Table III B-1, Specific QA/QC Requirements and Performance Standards for SW-846 Method 7470A and 7471A and are summarized below in Table III B-2 as "Required Analytical Deliverables (**YES**)". These routine reporting requirements should always be included as part of the laboratory deliverable for this method. It should be noted that although certain items are not specified as "Required Analytical Deliverables (**NO**)", these data are to be available for review during an audit and may also be requested on a client-specific basis.



Title: Quality Assurance and Quality Control Requirements and Performance Standards
for SW-846 Methods 7470A and 7471A, Mercury in Liquid Solid or Semi-Solid Waste

Table III B-2 Routine Analytical Reporting Requirements for SW-846 Methods 7470A/7471A

Parameter	Required Analytical Deliverable
Initial Calibration	NO
Initial Calibration Verification (ICV)	NO
Initial Calibration Blank (ICB)	NO
Low Level Calibration Check standard	NO
Continuing Calibration Verification (CCV)	NO
Continuing Calibration Blank (CCB)	NO
Method (Preparation) Blank	YES
Laboratory Control Standard (LCS)	YES
Field Matrix Spike Sample (MS)	YES, <u>only</u> if requested by the LSP
Field Matrix Duplicate (MD)	YES, <u>only</u> if requested by the LSP
Field Matrix Spike Duplicate (MSD)	YES, <u>only</u> if requested by the LSP
General Reporting	YES ¹

1. Non detected values must be reported with the sample-specific reporting limit.

4.0 Regulatory Limits for Mercury under 310 CMR 40.000

The most stringent (lowest) MCP Reportable Concentrations (RCs) and Method 1 Standards for mercury are as follows:

MCP Regulatory Criteria	Concentration
RC GW-1	0.001 mg/L (ppm)
RC S-1	20 mg/kg (ppm)
Groundwater Method 1 GW-3 Standard	1 µg/L (ppb)
Method 1 Soil Category S-1 & GW-1 Standard	20 µg/g (ppm)



Title: Sample Preservation, Container and Analytical Holding Time Specifications for Surface Water, Groundwater, Soil and Sediment Samples Analyzed for Metals by CVAA (SW-846 Methods 7470A/7471A)

Sample preservation, container and analytical holding time specifications for surface water, groundwater, soil, and sediment matrices for mercury analyses conducted in support of MCP decision-making are summarized below and presented in Appendix VII-A of CAM-VII A, Quality Assurance and Quality Control Guidelines for Sampling, Data Evaluation, and Reporting Activities for the Massachusetts Contingency Plan (MCP). Additional guidance may be found in SW-846, Chapter Three. Containers used for the collection of samples for the determination of Mercury should be pre-washed with detergent, acid washed and final rinsed with reagent grade water.

Sample Matrix	Container ¹	Preservation	Holding Time ²
Aqueous Total Mercury	500 mL glass or Polyethylene Bottle	HNO ₃ to pH < 2,	28 days
Aqueous Dissolved Mercury	500 mL glass or Polyethylene Bottle	Filter (0.45 µm) on site; or at the laboratory (prior to acid preservation) within 24 hours of collection HNO ₃ to pH <2,	28 days
Soil and Sediments	(1) 4-ounce glass jar with teflon-lined cap	Cool, 4°C	28 days
Concentrated Waste Samples	125 mL wide mouth glass or plastic	Cool to 4°C	28 days

1 The number of sampling containers specified is not a requirement. For specific analyses, the collection of multiple sample containers is encouraged to avoid resampling if sample is consumed or compromised.

2 From date of sample collection.



Title: MINIMUM FIELD QC SAMPLE FREQUENCY IN SUPPORT OF "PRESUMPTIVE CERTAINTY"
STATUS for Aqueous, Soil and Sediment Samples Analyzed for Mercury by CVAA
(SW-846 methods 7470A/7471A)

Parties seeking Presumptive Certainty for mercury analyses are required to provide the laboratory with the appropriate volumes and types of quality control (QC) samples as summarized in Table III B-1. Except for samples that require prior laboratory knowledge (i.e., Matrix Spikes, Matrix Spike Duplicates, drinking water samples, trip blanks/equipment blanks), field QC samples must be submitted to the laboratory as blind samples, using the containers and preservatives specified in CAM VII A, Appendix VII A-1, and Appendix III B-1 of this document.

MINIMUM FIELD QC SAMPLE FREQUENCY IN SUPPORT OF "PRESUMPTIVE CERTAINTY" STATUS

ANALYTES	METHOD (S)	QC ELEMENT	MEDIA		
			SOIL/SEDIMENT	AQUEOUS	DRINKING WATER
Mercury	SW-846 Method 7470/7471	Matrix Spike (MS)	1 per 20 samples ^d	1 per 20 samples ^d	1 per 20 samples ^d
		Field Duplicate	Not Mandatory ^a	Not Mandatory ^a	If analyte detected ^c
		Trip Blank	Not Mandatory ^a	Not Mandatory ^a	Not Mandatory ^a

^a Matrix Spikes must be selected that represent the most significant exposure points to human health and the environment. Separate MS required for sediment if both soil and sediment are sampled.

^b Samples designated "Not Mandatory" are not needed for Presumptive Certainty acceptance of analytical data. However, on a site and project-specific basis, the use of one or more of these and other QC elements (e.g. equipment rinsate blanks, etc.) may be advisable and/or necessary to demonstrate usability of the data, and/or to determine if the data are biased high due to contamination by sampling equipment/storage conditions.

^c Field Duplicate MUST be analyzed for any drinking water sample collected that mercury is detected in the primary sample above the reporting limit (RL). Contingency duplicate samples MUST be collected, consistent with the requirements of Appendix III B-1, for every trace metal sample, including drinking water samples for this purpose.